Invacare [®] Wheelchair Accessories			
Price List and Order Form Effective Octob		INVACARE	
For ease of ordering, contact Customer Service Toll Fre 800-333-6900 Fax: 800-678-4682	ee at:		8
To order Literature: 800-828-6282		www.invacare.com	Yes, you can.™
Company Name		Account #	
Phone #	P.O.#		Date

Ship To Name & Address __

Cyclical Lever Drive		Suggested List
Model#	Description	Price
LDR96LI	Cyclical Lever Drive - Left Hand For a 9000 Series Chair with 6" Casters	\$797.00
LDR98LI	Cyclical Lever Drive - Left Hand For a 9000 Series Chair with 8" Casters	797.00
LDR96RI	Cyclical Lever Drive - Right Hand For a 9000 Series Chair with 6" Casters	797.00
LDR98RI	Cyclical Lever Drive - Right Hand For a 9000 Series Chair with 8" Casters	797.00

Footnotes

- 1. Not available on a Super Low seat-to-floor height of 14 $1/2^{"}$.
- 2. Only available on 16",17",18" deep 9000 XDT.
- 3. Not available with 9 spoke rear wheels.
- 4. Not available with urethane tires.
- 5. Not available on a 9000 XT recliner.
- 6. Client weight not to exceed 250 lbs. on all chair models.
- 7. Seat-to-floor heights, see chart below for specific requirements.
- 8. Hemi seat-to-floor positions require the use of space-saver arms due to interference.

*Note: All specifications and dimensions are approximate.

Service and Maintenance Requirements

There are no service requirements for the CLD but the service requirements of the 9000 series wheelchair apply without change.

Specifications				
Weight Limit Mechanism Weight Added Width to Chair Stroke Length Grip Height Adjustment (from seat rail of 9000XT) Mechanism Color	250 lbs. 9.1 lbs. 1 ¹ /2" 15" - 16" 15" - 18" Black / Chrome			
Warranty	Two Years			
	Weight Limit Mechanism Weight Added Width to Chair Stroke Length Grip Height Adjustment (from seat rail of 9000XT) Mechanism Color	Mechanism Weight9.1 lbs.Added Width to Chair1 1/2"Stroke Length15" - 16"Grip Height Adjustment15" - 18"(from seat rail of 9000XT)Hechanism ColorMechanism ColorBlack / Chrome		

Installation

Right-handed or left-handed configuration. Dealer Installed Kit.

Caster	Seat To Floor Position	Caster Mtg	Rear Wheel Size Position	Rear Wheel	Space Saver Arm (9000 - order as required)
6"	15 ¹ /2"	top	20"	top	required
6"	16 ¹ /2"	middle	22"	top	required
6"	17 ¹ /2"	bottom	24"	top	required
6"	17 ¹ /2"	bottom	20"	bottom	no
8"	17 ¹ /2"	top	20"	bottom	no
8"	I7 ¹ /2"	top	24".	top	required
8"	18 ¹ /2".	middle	22"	bottom	no
8"	19 ¹ /2"	bottom	24"	bottom	no

Invacare Corporation w

on www.invacare.com



Yes, you can.[™]

USA One Invacare Way Elyria, OH USA 44035-4190 800-333-6900 Invacare is a registered trademark and Yes, you can. is a trademark of Invacare Corporation.

Specifications and prices are subject to change without notification. $\ensuremath{\textcircled{\sc 0}}$ 2001 Invacare Corporation